

	<b>Requirements for Downloading Software to Computer</b>	Document No	TD.01.DE.112
		Release Date	20.03.2025
		Revision Date	04.06.2025
		Revision No	01
		Page No:	1 / 2

- Minimum Hardware Requirements:
  - Minimum 4GB Ram
  - 1GHz CPU
  - Minimum i3 processor
  - 1366 x 768 Display Resolution
  - 1 USB Port
  - 1 GB free space
- Minimum Software Requirements:
  - Microsoft Windows8, Windows10, Windows11
  - .Net Framework 4.5 Full

The latest Fitting version below link:

[http://dosya.eartechnic.com/EarTuning2/EarTuning2\\_Setup.rar](http://dosya.eartechnic.com/EarTuning2/EarTuning2_Setup.rar)

[http://dosya.eartechnic.com/CorrectFit2/CorrectFit2\\_Setup.rar](http://dosya.eartechnic.com/CorrectFit2/CorrectFit2_Setup.rar)

**Note1:** Unzip the software file before setting up

**Note2:** Make sure the antivirus and firewall software are turned off during the download of EarTuning2 or CorrectFit2 Software.

**Note3:** Do not turn off or restart your computer during the download.

**Note4:**

- HiPro and HiPro2:
  - ✓ Make sure the original HiPro device and software used.
  - ✓ Configure HiPro

Prepared by	Controlled by	Approved by
Production Manager	Quality Manager	Quality Manager

	<b>Requirements for Downloading Software to Computer</b>	Document No	TD.01.DE.112
		Release Date	20.03.2025
		Revision Date	04.06.2025
		Revision No	01
		Page No:	2 / 2

### **User Roles and Permissions:**

The program allows users to have different levels of access. For example, a clinical expert can change all settings of the device, while a user can only adjust the volume. This is important to prevent unauthorized access and prevent accidental changes to the system.

### **Automatic Software Updates:**

Adaptation software must be maintained with continuous security patches and updates. Older software versions may be vulnerable to known vulnerabilities, so automatic updates or regular manual updates of the software should be provided.

### **Malware Detection:**

If users' devices or computers are infected with malware, the software can detect it by running a security scan and alert the user.

Prepared by	Controlled by	Approved by
Production Manager	Quality Manager	Quality Manager